


Summary

Name	AGRL4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9HBW9
Host	HEK293
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QKTTEFDTNS TDIALKVFFF DSYNMKHIHP HMNMDGDYIN IFPKRKAAYD SNGNVAFAFV YYKSIGPLLS SSDNFLKPKQ NYDNSEEEER VISSVISVSM SSNPPTLYEL EKITFTLSHR KVTDRYRSLC AFWNYSPDTM NGSWSSEGCE LTYSNETHTS CRCNHLTHFA ILMSSGPSIG IKDYNILTRI

image



Background

Alternative Names	Adhesion G protein-coupled receptor L4 (EGF, latrophilin and seven transmembrane domain-containing protein 1) (EGF-TM7-latrophilin-related protein) (ETL protein)
--------------------------	---

Note

Product Name: Recombinant Human ADGRL4 Protein
Catalog #: PHH6752



For Research Use Only , Not for Diagnostic Use