


Product Name: Recombinant Human FOLR3 Protein
Catalog #: PHH6758




Summary

Name	FOLR3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P41439
Host	HEK293
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QPRSARARTD LLNVCMNAKH HKTQPSPED ELYGQCSPWKK NACCTASTSQ ELHKDTSRLY NFNWDHCGKM EPTCKRHFIQ DSCLYECSPN LGPWIRQVNQ SWRKERILNV PLCKEDCERW WEDCRTSYTC KSNWHKGWNW TSGINECPAG ALCSTFESYF PTPAALCEGL WSHSFKVSNY SRGSGRCIQM WFDSAQGNPN EEVAKFYAAA MNAGAPSRGI IDS

image



SDS



WB

Background

Alternative Names	Folate receptor gamma (FR-gamma) (Folate receptor 3)
--------------------------	--

Note

Product Name: Recombinant Human FOLR3 Protein
Catalog #: PHH6758



For Research Use Only , Not for Diagnostic Use