


Product Name: Recombinant Human CADM4 Protein
Catalog #: PHH6767



Summary

Name	CADM4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NFZ8
Host	HEK293
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NVTVAEGGVA EITCRLHQYD GSIVVIQNPA RQTLFFNGTR ALKDERFQLE EFSPRRVRIR LSDARLEDEG GYFCQLYTED THHQIATLTV LVAPENPVVE VREQAVEGGE VELSCLVPRS RPAATLRWYR DRKELKGVSS SQENGKVWSV ASTVRFRVDR KDDGGIICE AQNQALPSGH SKQTQYVLVDV QYSPTARIHA SQAVVREGDT LVLTCAVTGN PRPNQIRWNR GNESLPERAE AVGETLTLPG LVSADNGTYT CEASNKHGHA RALYVLVVDYD PGAVVEAQT

image



Background

Alternative Names	Cell adhesion molecule 4 (Immunoglobulin superfamily member 4C) (IgSF4C) (Nectin-like protein 4) (NECL-4) (TSLC1-like protein 2)
--------------------------	--

Product Name: Recombinant Human CADM4 Protein
Catalog #: PHH6767



Note

For Research Use Only , Not for Diagnostic Use