


**Product Name: Recombinant Human CD69 Protein**  
**Catalog #: PHH6784**



## Summary

<b>Name</b>	CD69_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-hFc & C-His
<b>Accession #</b>	Q07108
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	44 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	QYTFSMPSDS HVSSCEDWV GYQRKCYFIS TVKRSWTSQAQ NACSEHGATL AVIDSEKDMN FLKRYAGREE HWVGLKKEPG HPWKWSNGKE FNNWFNVTGS DKCVFLKNTV VSSMECEKNL

## image



## Background

**Alternative Names** Early activation antigen CD69 (Activation inducer molecule) (AIM) (BL-AC/P26) (C-type lectin domain family 2 member C) (EA1) (Early T-cell activation antigen p60) (GP32/28) (Leukocyte surface antigen Leu-23) (MLR-3) (CD antigen CD69)

## Note

**Product Name: Recombinant Human CD69 Protein**  
**Catalog #: PHH6784**



---

For Research Use Only , Not for Diagnostic Use