


**Product Name: Recombinant Human CD7 Protein**  
**Catalog #: PHH6789**



## Summary

<b>Name</b>	CD7_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-His
<b>Accession #</b>	P09564
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	18 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	QSPHCTTVPV GASVNITCST SGGLRGIYLR QLGPQPQDII YYEDGVVPTT DRRFRGRIDF SGSQDNLIT MHLQLSDTG TYTCQAITEV NVYGSGLTLVL VTEEQSQGWH RCSDAPPRAS ALPAPPTGSA LPDPQTASAL PDPPAASALP

## image



## Background

**Alternative Names** T-cell antigen CD7 (GP40) (T-cell leukemia antigen) (T-cell surface antigen Leu-9) (TP41) (CD antigen CD7)

## Note

For Research Use Only , Not for Diagnostic Use