



Product Name: Recombinant Human ADAM8 Protein
Catalog #: PHH6807



Summary

Name	ADAM8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P78325
Host	HEK293
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CGNLFVERGE QCDCGPPEDC RNRCCNSTTC QLAEGAQCAH GTCCQECKVK PAGELCRPKK DMCDLEEFCD GRHPECPEDA FQENGTPCSG GYCYNGACPT LAQQCQAFWG PGGQAAEESC FSYDILPGCK ASRYRADMCG VLQCKGGQQP LGRAICIVDV CHALTTEGDT AYEPVPEGTR CGPEKVCWKG RCQDLHVYRS SNCSAQCHNH

image



Background

Alternative Names Disintegrin and metalloproteinase domain-containing protein 8 (ADAM 8) (EC 3.4.24.-) (Cell surface antigen MS2) (CD antigen CD156a)

Note

Product Name: Recombinant Human ADAM8 Protein
Catalog #: PHH6807



For Research Use Only , Not for Diagnostic Use