



Summary

Name	CHODL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	Q9H9P2
Host	HEK293
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CFADFKHPCY KMAYFHELSS RVSFQEARLA CESEGGVLLS LENEAEQKLI ESMLQNLTKP GTGISDGDFW IGLWRNGDGQ TSGACPDLYQ WSDGSNSQYR NWTDEPCG SEKCVVMYHQ PTANPGLGGP YLYQWNDRC NMKHNYICKY EPEINPTAPV EKPYLTNQPG DTHQNVVTE

image



Background

Alternative Names Chondrolectin (Transmembrane protein MT75)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CHODL Protein
Catalog #: PHH6809

