


Product Name: Recombinant Human PDP1 Protein
Catalog #: PHH6823



Summary

Name	PDP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	tag free
Accession #	Q9P0J1
Host	HEK293
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DSRAMLGVQE EDGSWSAVTL SNDHNAQNER ELERLKLEHP KSEAKSVVKQ DRLLGLLMPF RAFGDVKFKW SIDLQKRVIE SGPDQLNDNE YTKFIPPNYH TPPYLTAEPE VTYHRLRPQD

image



Background

Alternative Names	[Pyruvate dehydrogenase [acetyl-transferring]]-phosphatase 1, mitochondrial (PDP 1) (EC 3.1.3.43) (Protein phosphatase 2C) (Pyruvate dehydrogenase phosphatase catalytic subunit 1) (PDPC 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PDP1 Protein
Catalog #: PHH6823

