



Summary

Name	PEPA5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	P0DJD9
Host	HEK293
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IMYKVPLIRK KSLRRTLSEK GLLKDFLKKH NLNPARKYFP QWEAPTLVDE QPLENYLDME YFGTIGIGTP AQDFTVVFDK GSSNLWVPSV YCSSLACTNH NRFNPEDSST YQSTSETVSI TYGTGSMTGI LGYDTVQVGG ISDTNQIFGL SETEPGSFLY YAPFDGILGL AYPSISSGA TPVFDNIWNQ GLVSQDLFSV YLSADDKSGS VVIFGGIDSS YVTGSLNWWP VTVEGYWQIT VDSITMNGET IACAEGCQAI VDTGTSLLTG PTSPIANIQS DIGASENSDG DMVVSCSAIS SLPDIVFTIN GVQYPVPPSA YILQSEGSCI SGFQGMNVPT ESGELWILGD VFIRQYFTVF DRANNQVGLA PVA

image



Background

Alternative Names Pepsin A-5 (EC 3.4.23.1) (Pepsinogen-5)

Product Name: Recombinant Human PGA5 Protein
Catalog #: PHH6827



Note

For Research Use Only , Not for Diagnostic Use