



Product Name: Recombinant Human PSMB4 Protein
Catalog #: PHH6870



Summary

Name	PSB4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P28070
Host	HEK293
Species	Human
Molecular Mass	54 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ADGESFLGYV DMLGVAYEAP SLATGYGAYL AQPLLREVLE KQPVLSQTEA RDLVERCMRV LYYRDARSYN RFQIATVTEK GVEIEGPLST ETNWDIAHMI

image



Background

Alternative Names	Proteasome subunit beta type-4 (EC 3.4.25.1) (26 kDa prosomal protein) (HsBPROS26) (PROS-26) (Macropain beta chain) (Multicatalytic endopeptidase complex beta chain) (Proteasome beta chain) (Proteasome chain 3) (HsN3)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PSMB4 Protein
Catalog #: PHH6870

