


Product Name: Recombinant Human PPIC Protein
Catalog #: PHH6872



Summary

Name	PPIC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P45877
Host	HEK293
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KRGPSVTAKV FFDVRIGDKD VGRIVIGLFG KVPKTVENF VALATGEKGY GYKGSKFHRV IKDFMIQGGD ITTGDGTGGV SIYGETFPDE NFKLKHYGIG WVSMANAGPD TNGSQFFITL TKPTWLDGKH VVFGKVIDGM TVVHSIELQA TDGHDRPLTN CSIINSGKID VKTPFVVEIA

image



Background

Alternative Names Peptidyl-prolyl cis-trans isomerase C (PPIase C) (EC 5.2.1.8) (Cyclophilin C) (Rotamase C)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PPIC Protein
Catalog #: PHH6872

