



Product Name: Recombinant Human PTP4A2 Protein
Catalog #: PHH6942



Summary

Name	TP4A2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q12974
Host	HEK293
Species	Human
Molecular Mass	55 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HVLDWPFDDG APPPNQIVDD WLNLLKTKFR EEPGCCVAVH CVAGLGRAPV LVALALIECG MKYEDAVQFI RQKRRGAFNS KQLLYLEKYR PKMRLRFRDT

image



Background

Alternative Names Protein tyrosine phosphatase type IVA 2 (EC 3.1.3.48) (HU-PP-1) (OV-1) (PTP(CAAXII)) (Protein-tyrosine phosphatase 4a2) (Protein-tyrosine phosphatase of regenerating liver 2) (PRL-2)

Note

For Research Use Only , Not for Diagnostic Use