



Product Name: Recombinant Human VAMP7 Protein
Catalog #: PHH6980



Summary

Name	VAMP7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P51809
Host	HEK293
Species	Human
Molecular Mass	55 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GTTILAKHAW CGGNFLEVTE QILAKIPSEN NKLTYSHGNY LFHYICQDRI VYLCITDDDF ERSRAFNFLN EIKKRFQTTY GSRAQTALPY AMNSEFSSVL

image



Background

Alternative Names	Vesicle-associated membrane protein 7 (VAMP-7) (Synaptobrevin-like protein 1) (Tetanus-insensitive VAMP) (Ti-VAMP)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use