



Product Name: Recombinant Human TMPRSS11B Protein
Catalog #: PHH7001



Summary

Name	TM11B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86T26
Host	HEK293
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VEKTYYYQGD FHISGVTYND NCENAASQAS TNLSKDIETK MLNAFQNSSI YKEYVKSEVI KLLPNANGSN VQLQLKFKFP PAEGVSMRTK IKAKLHQMLK NNMASWNAV P ASIKLMEISK AASEMLTNNC CGRQVANSII TGNKIVNGKS SLEGAWPWQA SMQWKGRHYC GASLISSRWL LSAAHCFACK NNSKDWTVNF GIVVNKPYMT RKVQNIIFHE NYSSPGLHDD IALVQLAEEV SFTEYIRKIC LPEAKMKLSE NDNVVVTGWG TLYMNGSFPV ILQEDFLKII DNKICNASYA YSGFVTDTML CAGFMSGEAD ACQNDSGGPL AYPDSRNIWH LVGIVSWG DG CGKKNKPGVY TRVTSYRNWI

image



Background

Alternative Names Transmembrane protease serine 11B (EC 3.4.21.-) (Airway trypsin-like

Product Name: Recombinant Human TMPRSS11B Protein
Catalog #: PHH7001



protease 5)

Note

For Research Use Only , Not for Diagnostic Use