



Summary

| | |
|------------------------------|--|
| Name | ATS2_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-hFc & C-His |
| Accession # | O95450 |
| Host | HEK293 |
| Species | Human |
| Molecular Mass | 68 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | KKWSPCSKPC GGSQFTKYG CRRRLDHKMOV HRGFCAALSK PKAIRRACNP QECSQPVWVT GEWEPCSQTC GRTGMQVRSV RCIQPLHDNT TRSVHAKHCN DARPESRRAC SRELCPGRWR AGPWSQCSVT CGNGTQERPV LCRTADDSFG ICQEERPETA RTCRLGPCPR NISDPSKKSYP VVQWLSRPDP DSPIRKISSK GHCQGDKSIF CRMEVLSRYC SIPGYNKLCCK KSCNLYNNLT NVEGRIEPPP GKHNDIDVFM PTLPVPTVAM EVRPSPSTPL EVPLNASSTN ATEDHPETNA VDEPYKIHGL EDEVQPPNLI PRRPSPYEKT RNQRIQELID EMRKKEMLGK |

image



Background

Alternative Names

A disintegrin and metalloproteinase with thrombospondin motifs 2 (ADAM-TS 2) (ADAM-TS2) (ADAMTS-2) (EC 3.4.24.14) (Procollagen I N-proteinase)

Product Name: Recombinant Human ADAMTS2 Protein
Catalog #: PHH7018



(PC I-NP) (Procollagen I/II amino propeptide-processing enzyme)
(Procollagen N-endopeptidase) (pNPI)

Note

For Research Use Only , Not for Diagnostic Use