


Product Name: Recombinant Human MUC12 Protein
Catalog #: PHH7082



Summary

Name	MUC12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UKN1
Host	HEK293
Species	Human
Molecular Mass	11 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EESTASHSSQ DATGTIVLPA RSTTSVLLGE STTSPISSGS METTALPGST TTPGLSEKST TFHSSPRSPA TTLSPASTTS SGVSEESTTS

image



Background

Alternative Names Mucin-12 (MUC-12) (Mucin-11) (MUC-11)

Note

For Research Use Only , Not for Diagnostic Use