



Product Name: Recombinant Human SEMG1 Protein
Catalog #: PHH7112



Summary

Name	SEMG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P04279
Host	HEK293
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RLPSEFSQFP HGQKGQHYSG QKGKQQTESK GSFSIQYTYH VDANDHDQSR KSQQYDLNAL HKTTSQRHL GGSQQLLHNK QEGRDHDKSK GHFHRVVIHH KGGKAHRGTQ NPSQDQGNP SGKGISSQYS NTEERLWVHG LSKEQTSVSG AQKGRKQGGG QSSYVLQTEE LVANKQQRET KNSHQNKGHY QNVVEVREEH SSKVQTS LCP AHQDKLQHGS KDIFSTQDEL LVYNKNQHQT KNLNQDQQHG RKANKISYQS SSTEERRLHY GENGVQKDVS QSSIYSQTEE KAQGKSQKQI TIPSQEQEHS QKANKISYQS SSTEERRLHY GENGVQKDVS QRSIYSQTEK LVAGKSQIQA PNPKEPWHG

image



Background

Alternative Names	Semenogelin-1 (Cancer/testis antigen 103) (Semenogelin I) (SGI) [Cleaved
--------------------------	--

Product Name: Recombinant Human SEMG1 Protein
Catalog #: PHH7112



into: Alpha-inhibin-92; Alpha-inhibin-31; Seminal basic protein]

Note

For Research Use Only , Not for Diagnostic Use