


Product Name: Recombinant Human SMR3B Protein
Catalog #: PHH7116



Summary

Name	SMR3B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P02814
Host	HEK293
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QRGPRGPYPP GPLAPPQPFQ PGFVPPPPPP PYGPGRIPPP PPAPYGPQIF PPPPPQP

image



Background

Alternative Names Submaxillary gland androgen-regulated protein 3B (Proline-rich peptide P-B) (Proline-rich protein 3) [Cleaved into: Peptide P-A; Peptide D1A]

Note

For Research Use Only , Not for Diagnostic Use