


Product Name: Recombinant Human RNASE8 Protein
Catalog #: PHH7117



Summary

Name	RNAS8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TDE3
Host	HEK293
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KPKDMTSSQW FKTQHVQPSP QACNSAMSII NKYTERCKDL NTFLEPFSS VAITCQTPNI ACKNSCKNCH QSHGPMSTLM GELTSGKYPN CRYKEKHLNT PYIVACDPPQ QGDPGYPLVP VHLDKVV

image



SDS



WB

Background

Alternative Names Ribonuclease 8 (RNase 8) (EC 3.1.27.-)

Note

For Research Use Only , Not for Diagnostic Use