



Product Name: Recombinant Human MUC17 protein
Catalog #: PHH7135



Summary

Name	MUC17_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q685J3
Host	HEK293
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TQVTTSTEAS SSPPTAEVNS MPTSTPSEGS TPLTSMVSST MPVASSEAST LSTTPVDTST PVTTSSEASS SSTTPEGTSI PTSTPSEGST PLTNMPVSTR LVVSSEASTT STTPADSNTF VTTTSSEASSS STTAEGTSMP

image



Background

Alternative Names Mucin-17 (MUC-17) (Small intestinal mucin-3) (MUC-3)

Note

For Research Use Only , Not for Diagnostic Use