



Summary

Name	CHIA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BZP6
Host	HEK293
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YQLTCYFTNW AQYRPLGRF MPDNIDPCLC THLIYAFAGR QNNEITIEW NDVTLYQAFN GLKNKNSQLK TLLAIGGWNF GTAPFTAMVS TPENRQTFIT SVIKFLRQYE FDGLDFDWEY PGSRGSPQD KHLFTVLVQE MREAFEQEAK QINKPRLMVT AAVAAGISNI QSGYEIPQLS QYLDYIHVMT YDLHGSWEGY TGENSPYKY PTDTGSNAYL NVDYVMNYWK DNGAPAekli VGFPTYGHNF ILSNPSNTGI GAPTSGAGPA GPYAKESGIW AYEICTFLK NGATQGWDAP QEVPIYQGN VVVGYDNIKS FDIKAQWLKH NKFGGAMVWA IDLDDFTGTF CNQGFPLIS TLKKALGLQS ASCTAPAQPI EPITAAPSGS GNGSGSSSSG GSSGGSGFCA VRANGLYPVA NNRNAFWHCV NGVTYQQNCQ AGLVFDTS CCNWA

image



Background

Product Name: Recombinant Human CHIA Protein
Catalog #: PHH7141



Alternative Names

Acidic mammalian chitinase (AMCase) (EC 3.2.1.14) (Lung-specific protein TSA1902)

Note

For Research Use Only , Not for Diagnostic Use