


Product Name: Recombinant Human CYBB Protein
Catalog #: PHH7144



Summary

Name	CY24B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P04839
Host	HEK293
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FAVHHDEEKD VITGLKQKTL YGRPNWDNEF KTIASQHPNT RIGVFLCGPE ALAETLSKQS ISNSES GPRG

image



Background

Alternative Names	Cytochrome b-245 heavy chain (EC 1.-.-.-) (CGD91-phox) (Cytochrome b(558) subunit beta) (Cytochrome b558 subunit beta) (Heme-binding membrane glycoprotein gp91phox) (NADPH oxidase 2) (Neutrophil cytochrome b 91 kDa polypeptide) (Superoxide-generating NADPH oxidase heavy chain subunit) (gp91-1) (gp91-phox) (p22 phagocyte B-cytochrome)
--------------------------	---

Note

Product Name: Recombinant Human CYBB Protein
Catalog #: PHH7144



For Research Use Only , Not for Diagnostic Use