


Product Name: Recombinant Human SAA1 Protein
Catalog #: PHH7153



Summary

| | |
|------------------------------|--|
| Name | SAA1_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | N-MBP & C-His |
| Accession # | P0DJ18 |
| Host | HEK293 |
| Species | Human |
| Molecular Mass | 54 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | RSFFSFLGEA FDGARDMWRA YSDMREANYI GSDKYFHARG NYDAAKRGPG GVWAAEAISD ARENIQRFFG HGAEDSLADQ AANEWGRSGK DPNHFRPAGL PEKY |

image



Background

| | |
|--------------------------|---|
| Alternative Names | Serum amyloid A-1 protein (SAA) [Cleaved into: Amyloid protein A (Amyloid fibril protein AA); Serum amyloid protein A(2-104); Serum amyloid protein A(3-104); Serum amyloid protein A(2-103); Serum amyloid protein A(2-102); Serum amyloid protein A(4-101)] |
|--------------------------|---|

Note

Product Name: Recombinant Human SAA1 Protein
Catalog #: PHH7153



For Research Use Only , Not for Diagnostic Use