

Product Name: Recombinant Human SLC22A5 Protein
Catalog #: PHH7183



Summary

Name	S22A5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	O76082
Host	HEK293
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RVPDAANLSS AWRNHTVPLR LRDGREVPHS CRRYRLATIA NFSALGLEPG RDVDLGQLEQ ESCLDGWEFS QDVYLSTIVT EWNLVCEDDW

image



Background

Alternative Names Solute carrier family 22 member 5 (High-affinity sodium-dependent carnitine cotransporter) (Organic cation/carnitine transporter 2)

Note

For Research Use Only , Not for Diagnostic Use