


Product Name: Recombinant Human HTR3B Protein
Catalog #: PHH7187



Summary

Name	5HT3B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	O95264
Host	HEK293
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NEFVDIERYP DLPYVYVNSS GTIENYKPIQ VVSACSLETY AFPFDVQNCS LTFKSILHTV EDVDLAFLRS PEDIQHDKKA FLNDSEWELL

image



Background

Alternative Names 5-hydroxytryptamine receptor 3B (5-HT3-B) (5-HT3B) (Serotonin receptor 3B)

Note

For Research Use Only , Not for Diagnostic Use