


Product Name: Recombinant Human SLC6A4 Protein
Catalog #: PHH7240



Summary

Name	SC6A4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P31645
Host	HEK293
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QLSACEDGED CQENGVLQKV VPTPGDKVES GQISNGYS AV PSPGAGDDTR HSIPATTTTL VAE LHQGERE

image



Background

Alternative Names	Sodium-dependent serotonin transporter (SERT) (5HT transporter) (5HTT) (Solute carrier family 6 member 4)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use