


Product Name: Recombinant Human JAK3 Protein
Catalog #: PHH7269



Summary

Name	JAK3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P52333
Host	HEK293
Species	Human
Molecular Mass	59 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GLWNGAQLYA CQDPTIFEER HLKYISQLGK GNFGSVELCR YDPLGDNTGA LVAVKQLQHS GPDQQRDFQR EIQLKALHS DFIVKYRGVS YGPGRQSLRL VMEYLPSPGCL RDFLQRHRAR LDASRLLLYS SQICKGMEYL GSRRCVHRDL

image



Background

Alternative Names Tyrosine-protein kinase JAK3 (EC 2.7.10.2) (Janus kinase 3) (JAK-3) (Leukocyte janus kinase) (L-JAK)

Note

For Research Use Only , Not for Diagnostic Use