



Product Name: Recombinant Human RAC1 Protein
Catalog #: PHH7279



Summary

Name	RAC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P63000
Host	HEK293
Species	Human
Molecular Mass	62 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DGAVGKTCLL ISYTTNAFPG EYIPTVFDNY SANVMVDGKP VNLGLWDTAG QEDYDRLRPL SYPQTDVFLI CFSLVSPASF ENVRAKWYPE VRHHCNTPI ILVGTKLDLR DDKDTIEKLE EKKLTPITYP QGLAMAKEIG AVKYLECSAL TQRGLKTVFD

image



Background

Alternative Names Ras-related C3 botulinum toxin substrate 1 (EC 3.6.5.2) (Cell migration-inducing gene 5 protein) (Ras-like protein TC25) (p21-Rac1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RAC1 Protein
Catalog #: PHH7279

