



Product Name: Recombinant Human BANF1 Protein
Catalog #: PHH7412



Summary

Name	BAF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O75531
Host	HEK293
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MTTSQKHRDF VAEPMGKPV GSLAGIGEVL GKKLEERGF D KAYVVLGQFL VLKKDEDLFR EWLKDTCGAN AKQSRDCFCG LREWCD AFL

image



Background

Alternative Names Barrier-to-autointegration factor (Breakpoint cluster region protein 1)
[Cleaved into: Barrier-to-autointegration factor, N-terminally processed]

Note

For Research Use Only , Not for Diagnostic Use