



Summary

Name	EPB41_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P11171
Host	HEK293
Species	Human
Molecular Mass	57 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TVKAEVKKED EPPEQAEPEP TEAWKVEKTH IEVTVPTSNG DQTQKLAEKT EDLIRMRKKK RERLDGENIY IRHSNLMLED LDKSQEEIKK HHASISELKK NFMESVPEPR PSEWDKRLST

image



Background

Alternative Names	Protein 4.1 (P4.1) (4.1R) (Band 4.1) (EPB4.1) (Erythrocyte membrane protein band 4.1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use