


Product Name: Recombinant Human SPRR2E Protein
Catalog #: PHH7509



Summary

Name	SPR2E_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST & C-His
Accession #	P22531
Host	HEK293
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSYQQQCKQ PCQPPVCPT PKCEPCPPP KCPEPCPPK CPQPCPPQQC QQKCPPVTPS PPCQPKCPPK SK

image



Background

Alternative Names Small proline-rich protein 2E (SPR-2E) (Small proline-rich protein II) (SPR-II)

Note

For Research Use Only , Not for Diagnostic Use