


Product Name: Recombinant Human SP100 Protein
Catalog #: PHH7513



Summary

Name	SP100_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P23497
Host	HEK293
Species	Human
Molecular Mass	53 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KQSELPVTCG EVKGTLYKER FKQGTSKKCI QSEDKKWFTP REFEIEGDRG ASKNWKLSIR CGGYTLKVLV ENKFLPEPPS

image



Background

Alternative Names	Nuclear autoantigen Sp-100 (Nuclear dot-associated Sp100 protein) (Speckled 100 kDa)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use