



**Product Name: Recombinant Human MT3 Protein**  
**Catalog #: PHH7517**



## Summary

<b>Name</b>	MT3_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-His & N-SUMO
<b>Accession #</b>	P25713
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	18 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	GGCTCADSC KCEGCKCTSC KKSCCSCCPA ECEKCAKDCV CKGGEEAAEAE

## image



## Background

**Alternative Names** Metallothionein-3 (MT-3) (GIFB) (GIF) (Growth inhibitory factor) (Metallothionein-III) (MT-III)

## Note

For Research Use Only , Not for Diagnostic Use