



Product Name: Recombinant Human MT3 Protein
Catalog #: PHH7518



Summary

Name	MT3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST & C-His
Accession #	P25713
Host	HEK293
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDPETCPCPS GGSTCADSC KCEGCKCTSC KKSCCSCCPA ECEKCAKDCV CKGGEEAAEAE AEKCSCCQ

image



Background

Alternative Names	Metallothionein-3 (MT-3) (GIFB) (GIF) (Growth inhibitory factor) (Metallothionein-III) (MT-III)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use