


Summary

Name	CD151_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P48509
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ELKENLKDTM TKRYHQPGHE AVTSAVDQLQ QEFHCCGSNN SQDWRDSEWI RSQEAGGRVV PDSCCKTVVA LCGQRDHASN

image



Background

Alternative Names	CD151 antigen (GP27) (Membrane glycoprotein SFA-1) (Platelet-endothelial tetraspan antigen 3) (PETA-3) (Tetraspanin-24) (Tspan-24) (CD antigen CD151)
--------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------

Note

For Research Use Only , Not for Diagnostic Use