



**Product Name: Recombinant Human OASL Protein**  
**Catalog #: PHH7717**



## Summary

<b>Name</b>	OASL_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-MBP & C-His
<b>Accession #</b>	Q15646
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	69 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	FHSFQEAAKH HKDVLRLIWK TMWQSQDLLD LGLEDLRMEQ RVPDALVFTI QTRGTAEPIT VTIVPAYRAL GPLSPNSQPP PEVYVSLIKA CGGPGNFCPS FSELQRNFVK HRPTKLKSL RLVKHWYQQY VKARSPRANL PPLYALELLT IYAWEMGTEE DENFMLDEGF TTVMDLLLEY EVICIYWTKY YTLHNAIIED CVRKQLKKER PIILDPADPT

## image



## Background

<b>Alternative Names</b>	2'-5'-oligoadenylate synthase-like protein (2'-5'-OAS-related protein) (2'-5'-OAS-RP) (59 kDa 2'-5'-oligoadenylate synthase-like protein) (Thyroid receptor-interacting protein 14) (TR-interacting protein 14) (TRIP-14) (p59 OASL) (p59OASL)
--------------------------	--

**Product Name: Recombinant Human OASL Protein**  
**Catalog #: PHH7717**

---



**Note**

For Research Use Only , Not for Diagnostic Use