



Product Name: Recombinant Human LCE3C Protein
Catalog #: PHH7747



Summary

Name	LCE3C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST & C-His
Accession #	Q5T5A8
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QPPPSCPSPK CPPKSPAQCL PPPSSDCALS SGGCGPSSES GCCLSHHRHF RSHQCRRQRS NSCDRGSGQQ GGGSCRGHGS

image



Background

Alternative Names	Late cornified envelope protein 3C (Late envelope protein 15) (Small proline-rich-like epidermal differentiation complex protein 3A)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use