



Product Name: Recombinant Human LCE3B Protein
Catalog #: PHH7749



Summary

Name	LCE3B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST & C-His
Accession #	Q5TA77
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QPLPKCPSPK CPPKSSAQCL PPASSCCAPR PGCCGGPSSE GGCCLSHHRC CRSHRCRRQS SNSCDRGSGQ QDGASDCGYG

image



Background

Alternative Names Late cornified envelope protein 3B (Late envelope protein 14)

Note

For Research Use Only , Not for Diagnostic Use