


Product Name: Recombinant Human TBC1D13 Protein
Catalog #: PHH7965




Summary

Name	TBC13_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9NVG8
Host	HEK293
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PGIAKANMGV SREDVTFEDH PLNPNPDSRW NTYFKDNEVL LQIDKDVRRL CPDISFFQRA TDYPCLLLD PQNEFETLRK RVEQTTLKSQ TVARNRSGVT NMSSPHKNSV PSSLNEYEV

image



SDS



WB

Background

Alternative Names TBC1 domain family member 13

Note

For Research Use Only , Not for Diagnostic Use