


Product Name: Recombinant Human PEF1 Protein
Catalog #: PHH7989



Summary

| | |
|------------------------------|--|
| Name | PEF1_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | N-MBP & C-His |
| Accession # | Q9UBV8 |
| Host | HEK293 |
| Species | Human |
| Molecular Mass | 57 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | PGAAGQAPGA PPGSYYPGPP NSGGQYGSGL PPGGGYGGPA PGGPYGPPAG GGPYGHPNPG MFPSGTPGGP YGGAAPGGPY GQPPSSYGA QQPGLYGQGG APPNVDPEAY SWFQSVSDSH SGYISMKELK QALVNCNWSS |

image



Background

Alternative Names Peflin (PEF protein with a long N-terminal hydrophobic domain) (Penta-EF hand domain-containing protein 1)

Note

For Research Use Only , Not for Diagnostic Use