


Product Name: Recombinant Human IDH3A Protein
Catalog #: PHH8184



Summary

Name	IDH3A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P50213
Host	HEK293
Species	Human
Molecular Mass	56 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ANIMRMSDGL FLQKCREVAE SCKDIKFNEM YLDTVCLNMV QDPSQFDVLV MPNLYGDILS DLCAGLIGGL GVTPSGNIGA NGVAIFESVH GTAPDIAGKD MANPTALLLS AVMMMLRHMGL FDHAARIEAA

image



Background

Alternative Names	Isocitrate dehydrogenase [NAD] subunit alpha, mitochondrial (EC 1.1.1.41) (Isocitric dehydrogenase subunit alpha) (NAD(+)-specific ICDH subunit alpha)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use