


Product Name: Recombinant Human LDLRAD4 Protein
Catalog #: PHH8192



Summary

Name	LRAD4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15165
Host	HEK293
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NQSRRREDGL PQEGCLWPSD SAAPRLGASE IMHAPRSRDR FTAPSFQIRD RFSRFQPTY YVQHEIDLPP TISLSDGEEP PPYQGCTLQ LRDPEQQMEL NRESVRAPPN RTIFDSLID IAMYSGGPCP PSSNSGISAS TCSSNGRMEG PPPTYSEVMG HHPGASFLHH QRSNAHRGSR LQFQQNNAES TIVPIKGKDR

image



Background

Alternative Names Low-density lipoprotein receptor class A domain-containing protein 4

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human LDLRAD4 Protein
Catalog #: PHH8192

