



Product Name: Recombinant Human NMB Protein
Catalog #: PHH8220



Summary

Name	NMB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-9×His
Accession #	P08949-2
Host	HEK293
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LPEPRSRASK IRVHSRGNLW ATGHFMGKKS LEPSSPSPLG TAPHTSLRDQ RLQLSHDLLG ILLKKALGV SLRPPAPQIQ EAAGNTAEM TPIMGQTQQR GLDCAHPGKV LNTLLMAPS

image



Background

Alternative Names Neuromedin-B [Cleaved into: Neuromedin-B-32; Neuromedin-B]

Note

For Research Use Only , Not for Diagnostic Use