


Product Name: Recombinant Human PANK2 Protein
Catalog #: PHH8230



Summary

Name	PANK2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9BZ23
Host	HEK293
Species	Human
Molecular Mass	63 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VHWAAPPSLS SGLHRLFLR GTRIPSSTTL SPPRHDSL SL DGGTVNPPRV REPTGREAFG PSPASSDWLP ARWRNGRGGR PRARLCSGWT AAE EARRNPT LG LLLGRQRL LLRMGGGRLG APMERHGRAS ATSVSSAGEQ AAGDPEGRRQ EPLRRRASSA SVPVAVGASAE GTRRDRLGSY SGPTSVSRQR

image



Background

Alternative Names	Pantothenate kinase 2, mitochondrial (hPanK2) (EC 2.7.1.33) (Pantothenic acid kinase 2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PANK2 Protein
Catalog #: PHH8230

