


Product Name: Recombinant Human GNRH1 Protein
Catalog #: PHH8249



Summary

Name	GON1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His & N-SUMO
Accession #	P01148
Host	HEK293
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QHWSYGLRPG GKRDAENLID SFQEIVKEVG QLAETQRFEC TTHQPRSPLR DLKGALESLI EEETGQKKI

image



Background

Alternative Names	Progonadoliberin-1 (Progonadoliberin I) [Cleaved into: Gonadoliberin-1 (Gonadoliberin I) (Gonadorelin) (Gonadotropin-releasing hormone I) (GnRH-I) (Luliberin I) (Luteinizing hormone-releasing hormone I) (LH-RH I); GnRH-associated peptide 1 (GnRH-associated peptide I)]
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GNRH1 Protein
Catalog #: PHH8249

