


Product Name: Recombinant Human GNRH1 Protein
Catalog #: PHH8250



Summary

| | |
|------------------------------|---|
| Name | GON1_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-hFc & C-His |
| Accession # | P01148 |
| Host | HEK293 |
| Species | Human |
| Molecular Mass | 36 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | QHWSYGLRPG GKRDAENLID SFQEIVKEVG QLAETQRFEC TTHQPRSPLR DLKGALES LI EEETGQKKI |

image



Background

| | |
|--------------------------|--|
| Alternative Names | Progonadoliberin-1 (Progonadoliberin I) [Cleaved into: Gonadoliberin-1 (Gonadoliberin I) (Gonadorelin) (Gonadotropin-releasing hormone I) (GnRH-I) (Luliberin I) (Luteinizing hormone-releasing hormone I) (LH-RH I); GnRH-associated peptide 1 (GnRH-associated peptide I)] |
|--------------------------|--|

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GNRH1 Protein
Catalog #: PHH8250

