

Summary

Name	F234B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-9×His
Accession #	A2RU67
Host	HEK293
Species	Human
Molecular Mass	55 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLHSTWSRHL GSQGGDLSP LELADVNGDG LRDVLLSFVM SRNGSAVGVS RPAANLVCLS GMNGSTLWSS LLPEARDIT CLELMPGSLA ETICLVTGTH KMLSAFNATS GKAIWTLNPN YLSNGTLAAP VVVLPLDDED GVRDLVVLAI GELQPDLCFL LVSGRTGNPV GRPVKYNIVG VGNLIGPQVY ITTNGAVYIL FGFGNIQAVA LRDIFVQAQN RDSSPPLQI EEPEWEKRRS INLSELIDVY SDGVELLQMV KAPDSNCSNL LITTRQSLVL LRGQNLTYPYW ALRLQGLRSQ PTPGYFTDDQ TLDFLLQIQD GVGMMKMMVV DGDSGSIVWS YRAPCHMKET PATSAVTSQ KSVFLWAEG LSAASPNSDI ILGTEPPSLH HLYLLHPAFP SILLDLANTT GTVTASEVGI NDWLKDAFYV TRTTGPSSEG HPAALVSKL SLRWALMEGQ MAQLQESTPK IGRGELRRFL SRIKFVEAPY

image



Background

Product Name: Recombinant Human FAM234B Protein
Catalog #: PHH8254



Alternative Names

Protein FAM234B

Note

For Research Use Only , Not for Diagnostic Use