


Product Name: Recombinant Human RHOD protein
Catalog #: PHH8277



Summary

Name	RHOD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O00212
Host	HEK293
Species	Human
Molecular Mass	56 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AFPESYTPTV FERYMVNLQV KGKPVHLHIW DTAGQDDYDR LRPLFYPDAS VLLLCFDVTS PNSFDNIFNR WYPEVNHFCCK KVPIIVVGCK TDLRKDKSLV NKLRRNGLEP VTYHRGQEMA

image



Background

Alternative Names Rho-related GTP-binding protein RhoD (Rho-related protein HP1) (RhoHP1)

Note

For Research Use Only , Not for Diagnostic Use