



Product Name: Recombinant Human PDE6A Protein
Catalog #: PHH8284



Summary

Name	PDE6A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P16499
Host	HEK293
Species	Human
Molecular Mass	60 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AQNGLICNIM NAPAEDFFAF QKEPLDESGW MIKNVLSMPI VNKKEEIVGV ATFYNRKDGK PFDEMDETLM ESLTQFLGWS VLNPDTYESM NKLENRKDIF QDIVKYHVKC DNEEIQKILK TREVYGKEPW ECEEEELAEI LQAELPDADK

image



Background

Alternative Names	Rod cGMP-specific 3',5'-cyclic phosphodiesterase subunit alpha (GMP-PDE alpha) (EC 3.1.4.35) (PDE V-B1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use