


Product Name: Recombinant Human TMEM102 Protein
Catalog #: PHH8321



Summary

Name	TM102_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q8N9M5
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CGTPIREMWQ DCLGPPVPGA RDSIHRTSE ESKDWQSSV DQPHSYVTEH EAPVSLEKSP SDVSASESPQ

image



Background

Alternative Names	Transmembrane protein 102 (Common beta-chain associated protein) (CBAP)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use